



Recombinant Human Protein kinase C-binding protein NELL2 (NELL2)

Product Code	CSB-MP859098HU
Abbreviation	Recombinant Human NELL2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	Q99435
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥90% as determined by SDS-PAGE.
Sequence	LGVDPSLQIDVLTELELGESTTGVRQVPGLHNGTKAFLFQDTPRSIKASTATAE QFFQKLNRNKHEFTILVTLKQTHLNSGVILSIHHLDHRYLELESSGHRNEVRLHY RSGSHRPHTVEVPYILADDKWHKLSLAISASHLILHIDCNKIYERVVEKPSTDLP LGTTFWLGQRNNAHGYFKGIMQDVQLLVMPQGFIACPDNLNRTCPTCNDFHG LVQKIMELQDILAKTSAKLSRAEQRMNRLDQCYCERTCTMKGTTYREFESWID GCKNCTCLNGTIQCETLICPNPDCPLKSALAYVDGKCKECKSICQFQGRTYF EGERNTVYSSSGVVCVLYECKDQTMKLVESGCPALDCPESHQITLSHSCCKV CKGYDFCSEHRHNCMENSICRNLNDRAVCSCRDGFALREDNAYCEDIDECAE GRHYCRENTMCVNTPGSFMCICKTGYIRIDDYSCTEHDECITNQHNCDENALC FNTVGGHNCVCKPGYTGNNGTTCKAFCKDGCRRNGGACIAANVCACPQGFTGP SCETDIDECSDFVQCDSTRANCINLPGWYHCECRDGYHDNGMFSPSGESCE DIDECGTGRHSCANDTICFNLDGGYDCRCPHGKNCTGDCIHDGKVKHNGQIW VLENDRCVSCSCQNGFVMCRRMVDCENPTVDLFCCPECDPRLSSQCLHQN GETLYNSGDTWVQNCQQCRCLQGEVDCWPLPCPDVECEFSILPENECPPRC VTDPCQADTIRNDITKTCLDEMNVVRFTGSSWIKHGTCTLCQCKNGHICCSV DPQCLQEL
Research Area	Neuroscience
Source	Mammalian cell
Target Names	NELL2
Expression Region	22-816aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.



Tag Info	C-terminal 6xHis-tagged
Mol. Weight	91.1 kDa
Protein Length	Full Length of Mature Protein
Image	 (Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.
Typical	No
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.