




Recombinant Human Myosin-9 (MYH9), partial

Product Code	CSB-EP015303HUa0
Abbreviation	Recombinant Human MYH9 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	P35579
Form	Liquid or Lyophilized powder
Storage Buffer	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol.If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Alias	Cellular myosin heavy chain, type AMyosin heavy chain 9Myosin heavy chain, non-muscle IIaNon-muscle myosin heavy chain A ;NMMHC-ANon-muscle myosin heavy chain IIa ;NMMHC II-a ;NMMHC-IIA
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Sensitivity	Not Test
Purity	≥90% as determined by SDS-PAGE.
Sequence	AQQAADKYLYVDKNFINNPLAQADWAAKKLVWVPSDKSGFEPASLKEEVGEE AIVELVENGKKVKVNKDDIQKMNPCKFSKVEDMAELTCLNEASVLHNLKERY SGLIYTYSGLFVVPYKNLPYSEEIVEMYKGKKRHEMPPHIYAITDTAYRSM MQDREDQSILCTGESGAGKTENTKKVIQYLAYVASSHKSKKDQGELERQLLQA NPILEAFGNAKTVKNDNSSRFGKFIRI
Research Area	Cell Cycle
Source	E.coli
Target Names	MYH9
Expression Region	2-241aa
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	34.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.