





Recombinant Human Ephrin-B2 (EFNB2), partial (Active)

Product Code	CSB-MP007466HU1
Abbreviation	Recombinant Human EFNB2 protein, partial (Active)
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20?/-80?. The shelf life of lyophilized form is 12 months at -20?/-80?.
Uniprot No.	P52799
Form	Lyophilized powder
Storage Buffer	Lyophilized from a 0.2 µm sterile filtered PBS, 6% Trehalose, pH 7.4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Biological Activity	①Measured by its binding ability in a functional ELISA. Immobilized Nipah virus G (CSB-MP862323NDT1) at 2 µg/mL can bind Human EFNB2 . The EC50 is 0.6010-0.7267 ng/mL. ②Human EFNB2 captured on Protein A Chip can bind Recombinant Nipah virus G (CSB-MP862323NDT1) with an affinity constant of 0.358 nM as detected by MetaSPR Assay (WeSPRTM 200).
Purity	Greater than 95% as determined by SDS-PAGE.
Sequence	IVLEPIYWNSSNSKFLPGQGLVLYPQIGDKLDIICPKVDSKTVGQYEYYKVYMV DKDQADRCTIKKENTPLLNCAKPDQDIKFTIKFQEFSPNLWGLEFQKNKDYIIS TSNGSLEGLDNQEGGVCQTRAMKILMKVGQDASSAGSTRNKDPTRRPELEA GTNGRSSTTSPFVKPNPGSSTDGNSAGHSGNNILGSEVALFA
Source	Mammalian cell
Target Names	EFNB2
Expression Region	28-229aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4? for up to one week.
Tag Info	C-terminal hFc1-tagged
Mol. Weight	49.7 kDa
Protein Length	Partial
Image	




(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Activity

Measured by its binding ability in a functional ELISA. Immobilized Nipah virus G (CSB-MP862323NDT1) at 2 µg/ml can bind Human EFNB2. The EC₅₀ is 0.6010-0.7267 ng/mL.



Activity

Human EFNB2 captured on Protein A Chip can bind Recombinant Nipah virus G (CSB-MP862323NDT1) with an affinity constant of 0.358 nM as detected by MetaSPR Assay (WeSPRTM 200).

Endotoxin

Less than 1.0 EU/ug as determined by LAL method.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.