



# Recombinant Danio rerio RuvB-like 2 (ruvbl2)

<b>Product Code</b>	CSB-BP020599DIL
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P83571
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Danio rerio (Zebrafish) (Brachydanio rerio)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>MAAQVATTKV PEVRDITRIE RIGAHSHIRG LGLDDALEPR QVSQGMVGGQ  ASRRAGLIL EMIKDGQIAG RAVLIAGQPG TGKTAIAMGI AQLGPDTPF  TALAGSEIFS LEMSKTEALS QAFRKAIGVR IKEETEIIEG EVVEIQIDRP  ATGTGAKVVK LTLKTTEMET IYDLGTMIE SLSKERVQAG DVITIDKATG  KISKLGRSFT RARDYDAMGA QTQFVQCPEG ELQKRKEVVH TVSLHEIDVI  NSRTQGFLAL FSGDTGEIKS EVREQINAKV SEWREEGKAE IIPGVLFIDE  VHMLDIECFSLNRALESDDLSPVLIMATNR GITRIRGTNY QSPHGIPIDM  LDRLLIATT PYTEKETRQI LKIRCEEEDV ELSEEAHTVL TRIGQETSLR  YAIQLISTAG LVCRKRRGTE VQVEDIKRVY SLFLDEARSS QYMKEYQDSF  LFNETQTSQM DTS</p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	ruvbl2
<b>Protein Names</b>	Recommended name: RuvB-like 2 EC= 3.6.4.12 Alternative name(s): Reptin zReptin
<b>Expression Region</b>	1-463
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.