



Recombinant Human Matrix metalloproteinase-14 (MMP14), partial, Yeast

Product Code	CSB-YP014661HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50281
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Source	Yeast
Target Names	MMP14
Protein Names	Recommended name: Matrix metalloproteinase-14 Short name= MMP-14 EC= 3.4.24.80 Alternative name(s): MMP-X1 Membrane-type matrix metalloproteinase 1 Short name= MT-MMP 1 Short name= MTMMP1 Membrane-type-1 matrix metallopr
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal His-tagged Tag-Free
Protein Length	Extracellular domain
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.