



Recombinant ADP-ribose pyrophosphatase (nudF)

Product Code	CSB-MP016179EGX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P83842
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:H1 (strain CFT073 / ATCC 700928 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MLKPDNLPVT FGKNDVEIIA RETLYRGFFS LDLYRFRHRL FNGQMSHEVR REIFERGHAA VLLPFDPVRD EVVLIQIRI AAYDTSETPW LLEMVAGMIE EGESVEDVAR REAIEEAGLI VKRTPVLSF LSPGGTSER SSIMVGEVDA TTASGIHGLA DENEDIRVHV VSREQAYQWV EEGKIDNAAS VIALQWLQLH HQALKNEWA
Source	Mammalian cell
Target Names	nudF
Protein Names	Recommended name: ADP-ribose pyrophosphatase EC= 3.6.1.13 Alternative name(s): ADP-ribose diphosphatase ADP-ribose phosphohydrolase Short name= ASPPase Adenosine diphosphoribose pyrophosphatase Short name= ADPR-PPase
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.