



Recombinant RNA polymerase II-associated factor 1 homolog (C55A6.9)

Product Code	CSB-YP017382CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P90783
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPQESSDAKR IEPPrKVDFM LKPRFTNTVP DVPFDAKFMT CPFVPLGRFV EFQPAAIYRD YKHAVICDDD MGLNVDLIDL KKYDEDPIET EIDEKDNILL EDDGAALKIA KRSQQHSLV PWMRKTEYIS TEFNRFGVTA DRQETKLGYN LKKNNQQVEDM YRDKQSQIDA INKTFEDVRK PVKEHYSKKG VKAVEESFV PDFDHWKHLF AHVQFDGDTI TTEFEEEDER QQARESSVIK AMEFEDQKFA AVFVPTIGCL TSFMDDLELE RPFDEDMKYE FLLSREYTFK MEHLPPRDRD VFIMYHRNNV FQYNEVDCNV KMTRKPKMAL SRKSKLTLTY RNPSELEQKD MNKREAELYE QPKTRKQEIL EKIQEKKEEG GDSSDQSSDS DDDKPQKSRS DSSSDVSSDD DSPRKKEPTV DSDSD
Source	Yeast
Target Names	pafo-1
Protein Names	Recommended name: RNA polymerase II-associated factor 1 homolog
Expression Region	1-425
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.