



Recombinant Drosophila melanogaster Negative elongation factor E (Nelf-E)

Product Code	CSB-BP019521DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92204
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MVYIHFPNNL TEEEQMLQAK YQKLKKKKKA LQAHKAPKPE PESSLTLKRP TDARDAREVA RKLKSGAIP AIQKQTKQDQ TSFKRPKGQE RAKRSTSETT VASYQPFSST QNDVAQETII SEIKEEPRR QNLYQHFATE RDREERGMPE KVPMDTAQPE KPRAGNTIFV SGNKVTEDFL KKTFFNDYGTI VNVSMIEIEKS RGFVSFAKPE SADRAIAEIH GKNVNGINLQ VQLARRQPQI EPINDASSSA VWSSIAASKS QKGS HKDHRE MVQYDEDFLL
Source	Baculovirus
Target Names	Nelf-E
Protein Names	Recommended name: Negative elongation factor E
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.