



Recombinant Mesocricetus auratus V-type proton ATPase catalytic subunit A (ATP6V1A)

Product Code	CSB-MP002396MRG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P86205
Product Type	Recombinant Protein
Immunogen Species	Mesocricetus auratus (Golden hamster)
Purity	>85% (SDS-PAGE)
Sequence	VGHSELVGEI IRGVNVSALS RHKIMLPPRF SMVQVWPVRL PANHPLLTGQ RRTALVANTS NMPVAAREAS IYTGITLSEY FRWAEALREI SGRLAEMPAD SGYPAYLGAR KHFPSVNWLI SYSKALDEYY DKHFTEFVPL RTVGMLSNMI SFYDMARIKA DYAQLLEDMQ NAFR
Source	Mammalian cell
Target Names	ATP6V1A
Protein Names	Recommended name: V-type proton ATPase catalytic subunit A Short name= V-ATPase subunit A EC= 3.6.3.14 Alternative name(s): V-ATPase 69 kDa subunit Vacuolar proton pump subunit alpha
Expression Region	1-174
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.