



# Recombinant Splicing factor U2AF 65 kDa subunit (uaf-1)

<b>Product Code</b>	CSB-YP025408CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P90978
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSDHQDGMKL EDIERQFLDV AQREGGLEAI QPTTGPLENE ENLKSSTGGG GGEDDNRKK RKRSRSRDRD TRRRRSRDR GERRGGGGGG DRDRSRSRER RRGGGGRDEP RRRGGDDEAR SRREPEPQKP REPKKYRFWD VPPTGFETTT PMEYKNMQAA GQVPRGSVQS AVPVVGPSVT CQSRRLYVGN IPFGCNEEAM LDFFNQQMHL CGLAQAPGNP ILLCQINLDK NFAFIEFRSI DETTAGMAFD GINFMGQQLK VRRPRDYQPS QNTFDMNSRM PVSTIVVDSA NKIFIGGLPN YLTEDQVKEL LCSFGPLKAF SLNVDSQGNS KGYAFAEYLD PTLTDQAIAG LNGMQLGDKQ LVVQLACANQ QRHNTNLPNS ASAIAGIDLS QGAGRATEIL CLMNMVTEDE LKADDEYEEI LEDVRDECSK YGIVRSLEIP RPYEDHPVPG VGKVFVEFAS TSDCQRAQAA LTGRKFANRT VVTSYYDVK YHNRQF
<b>Source</b>	Yeast
<b>Target Names</b>	uaf-1
<b>Protein Names</b>	Recommended name: Splicing factor U2AF 65 kDa subunit Alternative name(s): U2 auxiliary factor 65 kDa subunit Short name= U2AF65 U2 snRNP auxiliary factor large subunit
<b>Expression Region</b>	1-496
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.