



Recombinant Drosophila melanogaster Mediator of RNA polymerase II transcription subunit 20 (MED20)

Product Code	CSB-EP013653DLU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91641
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MGVTILQPYP LPEGKSGAHI IDQLSKRLLA LGATHAGQFL VDCETFISTP QPHNGAPGRA VHVLNHSEYP ASTFSIIDNG TGKQVAIVAD NIFDLLMLKM TNTFTSKKQT KIESRGARFE YGDFVIKLGs VTMMEHFKGI LIEIEYKSCV ILAYCWEMIR EMLQGFLGIA VNKDFPSYFA PQTIMTAMGQ QQLHAKHNDI FEPMDTVKQY LEQFTNYRKH VTLMGGMGSG PGSQQVGPNV HMSPAVAGLH RS
Source	E.coli
Target Names	MED20
Protein Names	Recommended name: Mediator of RNA polymerase II transcription subunit 20 Alternative name(s): Mediator complex subunit 20 Protein trf-proximal dMED20 dTrap26 dTrfp
Expression Region	1-252
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.