



Recombinant Aspergillus flavus Uricase (uaZ)

Product Code	CSB-BP025648AVD
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q00511
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Aspergillus flavus
Purity	>85% (SDS-PAGE)
Sequence	SAVKAARYG KDNVRVYKVH KDEKTGVQTV YEMTVCVLE GEIETSYTKA DNSVIVATDS IKNTIYITAK QNPVTPPELF GSILGTHFIE KYNHIHAAHV NIVCHRWTRM DIDGKPHPHS FIRDSEEKRN VQVDVVEGKG IDIKSSLSGL TVLKSTNSQF WGFLRDEYTT LKETWDRILS TDVDATWQWK NFSGLQEVRS HVPKFDATWA TAREVTLKTF AEDNSASVQA TMYKMAEQIL ARQQLIETVE YSLPNKHIFE IDLSWHKGLQ NTGKNAEVFA PQSDPNGLIK CTVGRSSLKS KL
Source	Baculovirus
Target Names	uaZ
Protein Names	Recommended name: Uricase EC= 1.7.3.3Alternative name(s): Urate oxidase
Expression Region	2-302
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.