



Recombinant *Xenopus laevis* Tropomyosin alpha-1 chain (tpm1)

Product Code	CSB-YP024104XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q01173
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MDAIKKKMQM LKLDKENALD RAEQAEADKK GAEDKSKQLE DELVALQKKL KGTEDELDKY SEALKDAQEK LELSDKKATD AEGDVASLNR RIQLVEEELD RAQERLSTAL QKLEEAEEKAA DESERGMKVI ENRALKDEEK MELQEIQQLKE AKHIAEEADR KYEEVARKLV IIEGDLERAE ERAELSESKC AELEEEELKTV TNNLKSLEAQ AEKYSQKEDK YEEEIKVLTLD KLKEAETRAE FAERTVAKLE KSIDDLEDEL YAQKLKYKAI SEELDHALND MTSI
Source	Yeast
Target Names	tpm1
Protein Names	Recommended name: Tropomyosin alpha-1 chain Alternative name(s): Alpha-tropomyosin Tropomyosin-1
Expression Region	1-284
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.