



Recombinant *Saccharomyces cerevisiae* DNA-directed RNA polymerase III subunit RPC10 (RPC11)

Product Code	CSB-MP018352SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04307
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MLSFCPSCNN MLLITSGDSG VYTLACRSCP YEFPIEGIEI YDRKKLPRKE VDDVLGGGWD NVDQTKTQCP NYDTCGGESA YFFQLQIRSA DEPMTTFYKC VNCGHRWKEN
Source	Mammalian cell
Target Names	RPC11
Protein Names	Recommended name: DNA-directed RNA polymerase III subunit RPC10 Short name= RNA polymerase III subunit C10 Alternative name(s): DNA-directed RNA polymerases III 12.5 kDa polypeptide RNA polymerase III subunit C11
Expression Region	1-110
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.