



Recombinant *Saccharomyces cerevisiae* Peroxisome biogenesis factor 10 (PEX10)

Product Code	CSB-EP017794SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q05568
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MKNDNKLQKE ALMRLSQLRF PFADAPSIVQ AHQKDEQIQG LLIMKVTELC KLIKSQLFVN SYPKELSIFA KLLYLLFTTG RRGRTLGE EY VDLTYTNRKG TRLAGRLKMI VFAFAYPLCP YFITKLYKKI MKNNKESKIE DTESVAAFCK GLLDFILDVH MTLFYFKGAF YSISKRIFGM RYVFKHILSK NEANFREEGS QKYKVLGYIL LAQNVMKWYP VLTSTLGSWI YGRKRTNSI TRSSVGLQER SEHESIEGIP KESQLTHINL SDKNQLPFIP EASRKILCL MNMSDPSCAP CGHLFCWSCL MSWCKERPEC PLCRQHCQPQ EILVLRQ
Source	E.coli
Target Names	PEX10
Protein Names	Recommended name: Peroxisome biogenesis factor 10 Alternative name(s): Peroxin-10 Peroxisomal biogenesis factor 10 Peroxisome assembly protein 10
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.