



Recombinant Rabbit Acrosin (ACR)

Product Code	CSB-EP001183RB-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48038
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Purity	≥85% (SDS-PAGE)
Sequence	KDNA TCDGPCGLRF RQNPQGGFRV VGGQAAQQGA WPWMVSLQIF TPRNNRRYHA CGGVLLNAHW VLTAAHCFNN KQKVYEWRMV FGAQEIEYGT DKPVRPPLQE RYVEKVVTHD QNYNMTGEND IALLKITPPV PCGPFIFPGC LPNSKAGPPK AAQTCYVAGW GYVKENAPRP SPTLMEARVD LINLELCNST QWYNGRITAS NLCAGYPSGK IDTCQGDSGG PLMCRENQGE PFVVQGITSW GVGCAKAKRP GIYTATWPFL DWIASRIGSN ALRMIQPATP TPPTTRSPGV HQPPSAHPPW YFQHASGPPH PHPHPHPHPH PRPPQPAAQ APPPPPPPPP PPPPPPPPPP PPPPPPPPAS TKPPQALSFA KRLQQLVEVL KGGSTFSGAK RVDTA AAAAE ATEIPEVTLA S
Source	E.coli
Target Names	ACR
Protein Names	Recommended name: Acrosin EC= 3.4.21.10Cleaved into the following 2 chains: 1. Acrosin light chain 2. Acrosin heavy chain
Expression Region	17-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.