



Recombinant Malate dehydrogenase (mdh)

Product Code	CSB-MP013623CJN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0C890
Product Type	Recombinant Protein
Immunogen Species	Prosthecochloris vibrioformis (Chlorobium vibrioforme)
Purity	≥85% (SDS-PAGE)
Sequence	MKITVIGAGN VGATTAFRIA DKKLARELVL LDVVEGIPQG KGLDMYETGP VGLFDTKITG SNDYADTADS DIVIITAGLP RKPGMTREDL LMKNAGIVKE VTDNIMKHSK NPIIIIVVSNP LDIMTHVAWV RSGLPKERVI GMAGVLDAAR FRSFIAMELG VSMQDINACV LGGHGDAMVP VVKYTTVAGI PISDLLPAET IDKLVERTRN GGAEIVEHLK QGSAFYSPGS SVVEMVESIV LDRKRVLP CA VGLEGQYGID KTFVGVVPVKL GRNGVEQIYE INLDQADLDL LQKSAKIVDE NCKMLESTIG
Source	Mammalian cell
Target Names	mdh
Protein Names	Recommended name: Malate dehydrogenase EC= 1.1.1.37
Expression Region	1-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.