



Recombinant Cryptococcus neoformans var. neoformans serotype D Sorting nexin-3 (SNX3)

Product Code	CSB-BP022374CTL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0CR60
Product Type	Recombinant Protein
Immunogen Species	Cryptococcus neoformans var. neoformans serotype D (strain JEC21 / ATCC MYA-565) (Filobasidiella neoformans)
Purity	>85% (SDS-PAGE)
Sequence	MARMAGRPQT FDEMYSVPES FLEIEIRNPM THGIGRKMYT DYEIVCMTNI PAFKLRHSAV RRRYSDFEAF RDILERESTR VNIPPLPGKV RVFTNRFSD VIEQRREGLQ RFLEIVAGHP LLQTGSKVLC AFLQDPAWDK SQWI
Source	Baculovirus
Target Names	SNX3
Protein Names	Recommended name: Sorting nexin-3
Expression Region	1-144
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.