



Recombinant Pongo pygmaeus NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2 (NDUFA2)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-MP015623EXP |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P0CB80 |
| Product Type | Recombinant Protein |
| Immunogen Species | Pongo pygmaeus (Bornean orangutan) |
| Purity | >85% (SDS-PAGE) |
| Sequence | AAAAASRGI GAKLGLREIR IHLCQRSPGS QGVRDFIEKR YVELKKANPD LPILIRECSD VQPKLWARYA FGQEKNVPLN NFSADQVTRA LENVLKPEAS TED |
| Source | Mammalian cell |
| Target Names | NDUFA2 |
| Protein Names | Recommended name: NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 2 Alternative name(s): Complex I-B8 Short name= CI-B8 NADH-ubiquinone oxidoreductase B8 subunit |
| Expression Region | 2-103 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |