



Recombinant Human Serpin B10 (SERPINB10)

Product Code	CSB-MP021066HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P48595
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDSLATSINQ FALELSKKLA ESAQGKNIFF SWSISTSLT IVYLGAKGTT AAQMAQVLQF NRDQGVKCDP ESEKKRMEF NLSNSEEIHS DFQTLISEIL KPNDDYLLKT ANAIYGEKTY AFHNKYLEDM KTYFGAEPQP VNFVEASDQI RKDINSWVER QTEGKIQNLL PDDSVDSTTR MILVNALYFK GIWEHQFLVQ NTTEKPFRIN ETTSKPVQMM FMKKKLHIFH IEKPKAVGLQ LYYKSRDLSL LILLPEDING LEQLEKAITY EKLNEWTSAD MMELYEVQLH LPKFKLEDSY DLKSTLSSMG MSDAFSQSKA DFSGMSSARN LFLSNVFHKA FVEINEQGTE AAAGSGSEID IRIRVPSIEF NANHPFLFFI RHNKTNLILF YGRLCSP
Source	Mammalian cell
Target Names	SERPINB10
Protein Names	Recommended name: Serpin B10 Alternative name(s): Bomapin Peptidase inhibitor 10 Short name= PI-10
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.