



Recombinant Mouse Gastrin (Gast), partial

Product Code	CSB-EP009273MO-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P48757
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	SWKPRSLLQ DASSGPGTNE DLEQRQFNKL GSASHHRRQL GPQGPQHFLA DLSKKQRPRM EEEEEAYGWM DF
Source	E.coli
Target Names	Gast
Protein Names	Recommended name: GastrinCleaved into the following 3 chains: 1. Gastrin-71 Short name= 2. G71 3. Big gastrinAlternative name(s): Gastrin-34 Short name= G34 Gastrin
Expression Region	22-92
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.