



# Recombinant *Kluyveromyces lactis* Glucose-6-phosphate isomerase (RAG2), partial

<b>Product Code</b>	CSB-YP009717KBK
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P12341
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Kluyveromyces lactis</i> (strain ATCC 8585 / CBS 2359 / DSM 70799 / NBRC 1267 / NRRL Y-1140 / WM37) (Yeast) ( <i>Candida sphaerica</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Source</b>	Yeast
<b>Target Names</b>	RAG2
<b>Protein Names</b>	Recommended name: Glucose-6-phosphate isomerase Short name= GPI EC= 5.3.1.9 Alternative name(s): Phosphoglucose isomerase Short name= PGI Phosphohexose isomerase Short name= PHI
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Partial
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.