



Recombinant *Drosophila melanogaster* Proteasome subunit alpha type-1 (Pros35)

Product Code	CSB-EP018865DLU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P12881
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MFRNQYDSDV TVWSPQGRLH QVEYAMEAVK LGTATVGLKN KDYAVLVALC KPTSELSDTQ RKIIPIDDHL GISIAGLTAD ARVLSRYLRS ECLNYKHSYD TTYPVSRLIT NLGNKMQTTT QRYDRRPYGV GLLVAGYDER GPHIYQVTPS ATFFNCKANS IGSRSQSART YLEKLNKFL DSSKDEIIRH GIRAILGTLP TDEQKDGAGQ YDITVAIVGK DQPFTILSNK DSAKHVAIAK ENDNDTPRND DDDDRPSPPE EPAAGPRDPE VLVATEQRP
Source	<i>E.coli</i>
Target Names	Prosalph6
Protein Names	Recommended name: Proteasome subunit alpha type-1 EC= 3.4.25.1 Alternative name(s): PROS-Dm35 Proteasome 35 kDa subunit
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.