



Recombinant *Xenopus laevis* G2/mitotic-specific cyclin-B2 (ccnb2)

Product Code	CSB-BP004808XBE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P13351
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus laevis</i> (African clawed frog)
Purity	>85% (SDS-PAGE)
Sequence	MATRRAAIPR EADNILGGAM RSKVQMNSRR AALGEIGNKV TVRGKPPAVK QSSNAVAKPS KMAATKVANV KTKHVPVKPV VAEAAPKVPS PVPMDVSLKE EELCQAFSDA LTSVEDIDAD DGGNPQLCSD YVMDIYNYLK QLEVQQSVHP CYLEGKEINE RMRAILVDWL VQVHSRFQLL QETLYMGVAI MDRFLQVQPV SRSKLQLVGV TSLLIASKYE EMYTPEVADF VYITDNAYTA SQIREMEMII LRLLNFDLGR PLPLHFLRRA SKSCSADAEQ HTLAKYLMEL TLIDYEMVHI KPSEIAAAAL CLSQKILGQG TWGTTQHYYT GYTEGDLQLI MKHMAKNITK VNQNLTKHVA VRNKYASSKL MKISTLPQLM APLITELAAS LS
Source	Baculovirus
Target Names	ccnb2
Protein Names	Recommended name: G2/mitotic-specific cyclin-B2
Expression Region	1-392
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.