



# Recombinant Papio hamadryas Alcohol dehydrogenase 1B (ADH1B)

<b>Product Code</b>	CSB-EP001354ERJ-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P14139
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Papio hamadryas (Hamadryas baboon)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	STAGKVIKC KAAVLWEVKK PFSIEDVEVA PPKAYEVRİK MVAVGICRTD DHVVSIGNLVS PLPAILGHEA AGIVESVGEV VTTVKPGDKV IPLFTPQCGK CRVCKSPEGN YCVKNDLSNP RGTLQDGTTR FTCRGKPIHH FVGTSTFSQY TVVDENAVAK IDAASPLEKV CLIGCGFSTG YGSAVNVAKV TPGSVCAVFG LGGVGLSAVM GCKAAGAARI IAVDINKDKF AKAKELGATE CİNPQDYKKP IQEVLKEMTD GGVDVSFEVI GRLDTMMASL LCCHEACGTS VIVGVPPDSQ NLSINPMLLL TGRTWKGAVY GGFKSREGIP KLVADFMKK FSLDALITHV LPFEKINEGF DLLRSGKSIR TVLTF
<b>Source</b>	E.coli
<b>Target Names</b>	ADH1B
<b>Protein Names</b>	Recommended name: Alcohol dehydrogenase 1B EC= 1.1.1.1 Alternative name(s): Alcohol dehydrogenase subunit beta
<b>Expression Region</b>	2-375
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.