



Recombinant Chlorocebus aethiops Proto-oncogene serine/threonine-protein kinase mos (MOS)

Product Code	CSB-MP014733DSU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P10650
Product Type	Recombinant Protein
Immunogen Species	Chlorocebus aethiops (Green monkey) (Cercopithecus aethiops)
Purity	>85% (SDS-PAGE)
Sequence	MPSPLALRPY LPSEFSPSVD ARPCSSPSEL PAKLLLGATP PRAPRLPRRL AWCSIDWEQV CLLQRLGAGG FGSVYKATYH GVPVAIKQVN KCTKNRLASR RSFWAELNVA RLRHDNIVRV VAASTRTPAG SNSLGTIIME FGGNVTLHQV IYGAASHPEG DAGEPHCSTG GPLTLGKCLK YSLDVVNGLL FLHSQSIVHL DLKPANILIS EQDVCKISDF GCSEKLEDLL CFQTPLYPLG GTYTHRAPEL LKGEGVTPKA DIYSFAITLW QMTTKQAPYS GERQHILYAV VAYDLRPSLS AAVFQDSLPG QRLGDVIRRC WRPSAAQRPS ARPLLVDLTS LKAIEFG
Source	Mammalian cell
Target Names	MOS
Protein Names	Recommended name: Proto-oncogene serine/threonine-protein kinase mos EC= 2.7.11.1 Alternative name(s): Oocyte maturation factor mos Proto-oncogene c-Mos
Expression Region	1-346
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.