



Recombinant Human Fms-related tyrosine kinase 3 ligand (FLT3LG), partial

Product Code	CSB-BP008734HU1
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49771
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	TQDC SFQHSPISSD FAVKIRELSD YLLQDYPVTV ASNLQDEELC GGLWRLVLAQ RWMERLKTVA GSKMQGLLER VNTEIHFVTK CAFQPPPSCL RFVQTNISRL LQETSEQLVA LKPWITRQNF SRCLELQCQP DSSTLPPPWS PRPLEATAPT APQP
Source	Baculovirus
Target Names	FLT3LG
Protein Names	Recommended name: Fms-related tyrosine kinase 3 ligand Short name= Flt3 ligand Short name= Flt3L Alternative name(s): SL cytokine
Expression Region	27-184
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.