



Recombinant *Saccharomyces cerevisiae* Coproporphyrinogen-III oxidase (HEM13)

Product Code	CSB-EP005910SVG-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11353
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPAPQDPRNL PIRQQMEALI RRKQAEITQG LESIDTVKFH ADTWTRGNDG GGGTSMVIQD GTTFEKGGVN VSVVYQQLSP AAVSAMKADH KNLRLPEDPK TGLPVTDGVK FFACGLSMVI HPVNPHAPTT HLNRYRFETW NQDGTPQTWW FGGGADLTPS YLYEEDGQLF HQLHKDALDK HDTALYPRFK KWCDEYFYIT HRKETRGIGG IFFDDYDERD PQEILKMVED CFDAFLPSYL TIVKRRKDMP YTKEEQWQA IRRGRYVEFN LIYDRGTQFG LRTPGSRVES ILMSLPEHAS WLYNHHPAPG SREAKLLEVT TKPREWVK
Source	E.coli
Target Names	HEM13
Protein Names	Recommended name: Coproporphyrinogen-III oxidase Short name= COX Short name= Coprogen oxidase Short name= Coproporphyrinogenase EC= 1.3.3.3
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.