



# Recombinant Saccharopolyspora erythraea Acyl carrier protein (acpP)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-EP015636STB-B   |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P11830  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Saccharopolyspora erythraea (strain ATCC 11635 / DSM 40517 / JCM 4748 / NBRC 13426 / NCIMB 8594 / NRRL 2338)  |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | MDRKEIFERI EQVLAEQLGI PAEQITEEAD LREDLGMDSL DLVELVSALE<br>DEVGMRVEQS QLEGIETVGH VMELTLDLVA RLATASAADK PEAAS   |
| <b>Source</b>            | E.coli  |
| <b>Target Names</b>      | acpP  |
| <b>Protein Names</b>     | Recommended name: Acyl carrier protein Short name= ACP  |
| <b>Expression Region</b> | 1-95  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | full length protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.  |