



Recombinant Mouse Tumor necrosis factor receptor superfamily member 3 (Ltbr), partial

Product Code	CSB-BP013227MO
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P50284
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	QLVPPYRIENQTCWDQDKEYYEPMHVCCSRCPGGEFVFAVCSRSQDTVCK TCPHNSYNE HWNHLSTCQLCRPCDIVLGFEEVAPCTSDRKAECRCQPGMSCVYLDNECVH CEEERLVLC QPGTEAEVTDEIMDTDVNCVPCKPGHFQNTSSPRARCQPHTRCEIQGLVEAA PGTSYSYSDT ICKNPPEPGAMLL
Source	Baculovirus
Target Names	Ltbr
Protein Names	Recommended name: Tumor necrosis factor receptor superfamily member 3 Alternative name(s): Lymphotoxin-beta receptor
Expression Region	31-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Extracellular domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.