



# Recombinant Sheep Thyroxine-binding globulin (SERPINA7)

<b>Product Code</b>	CSB-BP021063SH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P50450
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Ovis aries (Sheep)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	APPN SCEGKITSCL SPQQNATLYK MSSINADFAF NLYRKVTVET PDQNIFFSPV SISAGLAMLS LGACSSSTQTQ ILESLGFNLT DTPMAEIQQG FQHLICSLNF PKKELELQMG NALFIGKQLK PLEKFLDDVK NLYETEVFST DFSNVSAQQ EINSHVERQT KGKIVGLIQD LKPNTITVLV NYLCFKAQWA NPFDPKTEE GSSFLVDKTT TVQVPMMHQV DQYYHLVDTE LNCTVLQMDY SKNALALFVL PKEGQMEGVE GAMSSKILKK WNRLLQKGWV NLFVPKFSIS ATYDLGGILL KMGIQDAFAD NADFSGLTKD NGLKVSNAVAH KAMFYIGEKG TEAIPVRFLL NQPETLLHP IIQFDRSFL LILEKNTRSI LFLGKVDPT EV
<b>Source</b>	Baculovirus
<b>Target Names</b>	SERPINA7
<b>Protein Names</b>	Recommended name: Thyroxine-binding globulin Alternative name(s): Serpin A7 T4-binding globulin
<b>Expression Region</b>	17-412
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.