



# Recombinant *Saccharomyces cerevisiae* Geranylgeranyl transferase type-2 subunit beta (BET2)

<b>Product Code</b>	CSB-BP019235SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P20133
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSGSLTLLKE KHIRYIESLD TKKHNFEYWL TEHLRLNGIY WGLTALCVLD SPETFVKEEV ISFVLSCWDD KYGAFAPFPR HDAHLLTTLS AVQILATYDA LDVLGKDRKV RLISFIRGNQ LEDGSFQGDR FGEVDTRFVY TALSALSILG ELTSEVVDPA VDFVLKCYNF DGGFGLCPNA ESHAAQAFTC LGALAIANKL DMLSDDQLEE IGWWLCERQL PEGGLNGRPS KLPDVCYSWW VLSSLAIGR LDWINEKLT EFILKCQDEK KGGISDRPEN EVDVFHTVFG VAGLSLMGYD NLVPIDPIYC MPKSVTSKFK KYPYK
<b>Source</b>	Baculovirus
<b>Target Names</b>	BET2
<b>Protein Names</b>	Recommended name: Geranylgeranyl transferase type-2 subunit beta EC= 2.5.1.60 Alternative name(s): Geranylgeranyl transferase type II subunit beta Short name= GGTase-II-beta Type II protein geranyl-geranyltransferase subunit beta
<b>Expression Region</b>	1-325
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.