



Recombinant Fundulus heteroclitus L-lactate dehydrogenase B chain (ldhb)

Product Code	CSB-MP012841FDK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20373
Product Type	Recombinant Protein
Immunogen Species	Fundulus heteroclitus (Killifish) (Mummichog)
Purity	>85% (SDS-PAGE)
Sequence	MSSVLQKLIT PLASSSAEPP RNKVTVVGVG QVGMACAVSI LLRDLCDELA LVDVMEDRLK GEMMDLQHGL LFLKTSKVVA DKDYAVTANS RLVVVTAGVR QQEGESRLNL VQRNVNVFKC IIPQIIKYSP NCTILVVSNP VDVLTYYVTWK LSGLPKHRVI GSGTNLDSAR FRYMMAERLG IHASAFNGWV LGEHGDTSPV VWSGANVAGV SLQKLNPEIG TDGDKEQWKA THKAVVDSAY EVIKLKGYTN WAIGFSVADL TESIVKNLSR VHPVSTMVKD MFGIGEEVFL SLPCVLNGSG VGSVVNMTLT DAEVAQLKKS ADTLWGIQKD LKDL
Source	Mammalian cell
Target Names	ldhb
Protein Names	Recommended name: L-lactate dehydrogenase B chain Short name= LDH-B EC= 1.1.1.27
Expression Region	1-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.