



Recombinant Enterobacteria phage T3 N-acetylmuramoyl-L-alanine amidase (3.5)

Product Code	CSB-EP017863EDY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20331
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage T3 (Bacteriophage T3)
Purity	>85% (SDS-PAGE)
Sequence	AKVQFKPRA TTEAIFVHCS ATKPSQNVGV REIRQWHKEQ GWLDVGYHFI IKRDGTVEAG RDELAVGSHA KGYNHNSIGV CLVGGIDDKG KFDANFTPAQ MQSLRSLLVT LLAKYEGSVL RAHHDVAPKA CPSFDLKRWW EKNELVTSR G
Source	E.coli
Target Names	3.5
Protein Names	Recommended name: N-acetylmuramoyl-L-alanine amidase EC= 3.5.1.28 Alternative name(s): T3 lysozyme
Expression Region	2-151
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.