



Recombinant Dictyostelium discoideum Ras-related protein Rab-8A (rab8A)

Product Code	CSB-BP019222DKK
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P20790
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MTSPATNKSA AYDYLIKLLL IGDSGVGKSC LLLRFSEDSF TPSFITTIGI DFKIRTIELE GKRIKLIQWD TAGQERFRTI TTAYYRGAMG ILLVYDVTDE KSGFNIRNWI RNIEQHATDS VNKMLIGNKC DMAEKKVVDS SRGKSLADEY GIKFLETSK NSINVEEAFI SLAKDIKKRM IDTPNEQPQV VQPGTNLGN NNKKKACC
Source	Baculovirus
Target Names	rab8A
Protein Names	Recommended name: Ras-related protein Rab-8A Alternative name(s): GTP-binding protein SAS1
Expression Region	1-208
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.