



# Recombinant Daucus carota Tubulin beta-1 chain (TUBB1)

<b>Product Code</b>	CSB-EP025319DIR
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P20364
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Daucus carota (Wild carrot)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	FQVCHSLGGG TGSGMGTLII SKIREEYPDR MMLTFSVFPS PKVSDTVVEP YNATLSVHQL VENADVCMVL DNEALYDICF RTLKLTTPSF GDLNHLISAT MSGVTCSLRF PGQLNSDLRK LAVILIPFPR LHFFMVGFAP LTSRGSQQYR SLSVPELTQQ MWDSKNMMCA ADPRHGRYLT ASAMFREKMS TKDLDEQMIN VQKNSSYFV EWIPNNVKST VCDIPPTGLK MASTFIGNST SIQEMFRRVS EQFTAMFRRK AFLHWYTGEG MDEMEFTEAE SNMNDLVSEY QQYQDATADE EEYYEDEEEE EAQGM
<b>Source</b>	E.coli
<b>Target Names</b>	TUBB1
<b>Protein Names</b>	Recommended name: Tubulin beta-1 chain Alternative name(s): Beta-1-tubulin
<b>Expression Region</b>	1-315
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.