



# Recombinant *Notophthalmus viridescens* Retinoic acid receptor alpha (RARA)

<b>Product Code</b>	CSB-EP019338NEO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P18514
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Notophthalmus viridescens</i> (Eastern newt) ( <i>Triturus viridescens</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASNGGSCPS SGGHMNGYPV PHYAFFFFPHM LGGLSPPGSL AGIPHPLPVS AYSTPSPATI ETQSTSSEEI VSPSPSPPL PRIYKPCFVC QDKSSGYHYG VSACEGCKGF FRSIQKNMV YTCHRDKTCI INKVTRNRCQ YCRLQKCFEV GMSKESVRND RNKKKKQEAP KQECTESYII TPEVEDLVEK VRKAHQETFP ALCQLGKYTT NNSSEERVSL DIDLWDFSE LSTKCIKTV EFAKQLPGFT TLTIADQITL LKAACLDILI LRICTRYTPD QDTMTFSDGL TLNRTQMHNA GFGPLTDLVF AFANQLLPLE MDDAETGLLS AICLICGDRQ DLEQPDKVDK LQEPILLEALK IYVRKRRPNK PHMFPKMLMK ITDLRSISAK GAERVITLKM EIPGSMPLI QEMLENSEGL DSLTGQPPRA SSLAPPPGSC SPSLSPSSNR SSPTSHSP
<b>Source</b>	E.coli
<b>Target Names</b>	RARA
<b>Protein Names</b>	Recommended name: Retinoic acid receptor alpha Short name= RAR-alpha Alternative name(s): Nuclear receptor subfamily 1 group B member 1
<b>Expression Region</b>	1-458
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.