



Recombinant Anas platyrhynchos Alpha-enolase (ENO1)

Product Code	CSB-MP007670BZD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P19140
Product Type	Recombinant Protein
Immunogen Species	Anas platyrhynchos (Mallard) (Anas boschas)
Purity	>85% (SDS-PAGE)
Sequence	SILKIHARE IFDSRGNPTV EVDLYTNKGL FRAAVPSGAS TGIYEALRLR DNDKTRYMGK GVSKAVEHIN KTIAPALISK NVNVVEQDKI DKLMLDMDGS ENKSKFGANA ILGVSLAVCK AGAAEKGVPL YRHIADLAGN PEVILPVPAF NVINGGSHAG NKLAMQEFMI PPCGADSFKE AMRIGAEVYH NLKNVIKEKY GKDATNVGDE GGFAPNILEN KEALELLKTA IGKAGYSDKV VIGMDVAASE FYRDGKYDLD FKSPDDPSRY ISPDQLADLY KGFVKNYPVV SIEDPFDQDD WGAWKKFTGS VGIQVVGDDL TVTNPKRIAK AVEEKACNCL LLKVNQIGSV TESLQACKLA QSNWGWVMVS HRSGETEDTF IADLVVGLCT GQIKTGAPCR SERLAKYNQL LRIIEELGSK ARFAGRNFNRN PRIN
Source	Mammalian cell
Target Names	ENO1
Protein Names	Recommended name: Alpha-enolase EC= 4.2.1.11 Alternative name(s): 2-phospho-D-glycerate hydro-lyase Tau-crystallin
Expression Region	2-434
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.