




# Recombinant Rat Plasminogen activator inhibitor 1 (Serpine1)

<b>Product Code</b>	CSB-EP021081RA
<b>Relevance</b>	Serine protease inhibitor. This inhibitor acts as 'bait' for tissue plasminogen activator, urokinase, protein C and matriptase-3/TMPRSS7. Its rapid interaction with PLAT may function as a major control point in the regulation of fibrinolysis .
<b>Abbreviation</b>	Recombinant Rat Serpine1 protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P20961
<b>Alias</b>	Endothelial plasminogen activator inhibitor Serpin E1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Purity</b>	≥ 90% as determined by SDS-PAGE.
<b>Sequence</b>	SPLPESHTAQQATNFGVKVFQHVQASKDRNVVFSYGVSSVLAMLQLTTAG KTRQQIQDAMGFNISERGTAPALRKLSKELMGSWNKNEISTADAIFVQRDLELV QGFMPHFFKLFRTTVKQVDFSEVERARFIINDWVERHTKGMISDLLAKGAVNE LTRLVLVNALYFNGQWKTPFLEASTHQRLFHKSDGSTISVPMMAQNNKFNYTE FTTPDGHEYDILELPYHGETLSMFIAAPFEKDVPLSAITNILD AELIRQWKSNT RLPRLLILPKFSLETEVDLRGPLEKLGMTDIFSSTQADFTSLSDQEQLSVAQAL QKVKIEVNESGTVASSSTAILVSARMAPTEMVLD R SFLVVRHNPTETILFMGQ LMEP
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	Serpine1
<b>Protein Names</b>	Recommended name: Plasminogen activator inhibitor 1 Short name= PAI Short name= PAI-1 Alternative name(s): Endothelial plasminogen activator inhibitor Serpin E1
<b>Expression Region</b>	24-402aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-SUMO-tagged
<b>Mol. Weight</b>	58.6kDa
<b>Protein Length</b>	Full Length of Mature Protein
<b>Image</b>	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

### Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.