



# Recombinant Human 40S ribosomal protein S9 (RPS9)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP020482HU  |
| <b>Storage</b>           | Store at -20°C, for extended storage, conserve at -20°C or -80°C.   |
| <b>Uniprot No.</b>       | P46781  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Homo sapiens (Human)  |
| <b>Purity</b>            | ≥85% (SDS-PAGE)   |
| <b>Sequence</b>          | PVARSWVCR KTYVTPRRPF EKSRLDQELK LIGEYGLRNK REVWRVKFTL<br>AKIRKAAREL LTLDEKDP RR LFEGNALLRR LVRIGVLDEG KMKLDYILGL<br>KIEDFLERRL QTQVFKLGLA KSIHHARVLI RQRHIRVRKQ VVNIPSFIVR<br>LDSQKHIDFS LRSPYGGGRP GRVKRKNACK GQGGAGAGDD EEED  |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | RPS9  |
| <b>Protein Names</b>     | Recommended name: 40S ribosomal protein S9  |
| <b>Expression Region</b> | 2-194   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Target Details</b>    | Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit. Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins. This gene encodes a ribosomal protein that is a component of the 40S subunit. The protein belongs to the S4P family of ribosomal proteins. It is located in the cytoplasm. Variable expression of this gene in colorectal cancers compared to adjacent normal tissues has been observed, although no correlation between the level of expression and the severity of the disease has been found. As is typical for genes encoding ribosomal proteins, multiple processed pseudogenes derived from this gene are dispersed through the genome. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.   |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.