



Recombinant Avena sativa Tubulin beta-1 chain (TUBB1)

Product Code	CSB-YP025319DPO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P25862
Product Type	Recombinant Protein
Immunogen Species	Avena sativa (Oat)
Purity	>85% (SDS-PAGE)
Sequence	AVLMDLEPGT MDSVRTGPYG QIFRPDNFVF GQSGAGNNWA KGHYTEGAEL IDSVLDVVRK EAENCDCLQG FQVCHSLGGG TGSGMGTHLLI SKIREEYPDR MMLTFSVFPS PKVSDTVVEP YNATLSVHQL VENADECMLV DNEALYDICF RTLKLTTPSF GDLNHLISAT MSGVTCCLRF PGQLNSDLRK LAVNLIPFPR LHFFMVGAFAP LTRGSQQR ALTVPCLTQQ MWDSKNMMCA ADPRHGRLT ASAMFRGKMS TKEVDEQMIN VQKNSSYFV EWIPNNVKSS VCDIPPTGLS MASTFIGNST SIQEMFRRVS EQFTAMFRRK AFLHWYTGEG MDEMEFTEAE SNMNDLVSEY QQYQDATADE EGEYEDEEED LQAEDM
Source	Yeast
Target Names	TUBB1
Protein Names	Recommended name: Tubulin beta-1 chain Alternative name(s): Beta-1-tubulin
Expression Region	1-386
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.