



Recombinant Human Beta-crystallin B3 (CRYBB3)

Product Code	CSB-MP006015HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26998
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAEQHGAPEQ AAAGKSHGDL GGSYKVILYE LENFQGKRCE LSAECPSLTD SLLEKVGSIQ VESGPWLAFE SRAFRGEQFV LEKGDYPRWD AWSNSRSDSDS LLSLRPLNID SPHHLHLFE NPAFSGRKME IVDDDVPSLW AHGFQDRVAS VRAINGTWVG YEFPGYRGRQ YVFERGEYRH WNEWASQPQ LQSVRRIRDQ KWHKRGRFPS S
Source	Mammalian cell
Target Names	CRYBB3
Protein Names	Recommended name: Beta-crystallin B3 Alternative name(s): Beta-B3 crystallin
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.