



Recombinant Pisum sativum Tubulin beta-3 chain (TUBB3)

Product Code	CSB-MP025323EWE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29502
Product Type	Recombinant Protein
Immunogen Species	Pisum sativum (Garden pea)
Purity	>85% (SDS-PAGE)
Sequence	GQCGNQIGSK FWEVVCDEHG IDPTGRYVGN SDLQLERVNV YYNEASCGRF VPRAILMDLE PGTMDSVRTG PYGQIFRPDN FVFGQSGAGN NWAKGHYTEG AELIDSVLDV VRKEAENCDC LQGFQVCHSL GGGTGSGMGMT LLISKIREEY PDRMMLTFSV FPSPKVSDTV VEPYNATLSV HQLVENADEC MVL DNEALYD ICFRTLKLT T PSFGDLNHLI SATMSGVTCC LRFPGQLNSD LRKLAVNLIP FPR LHFFMVG FAPLTSRGSQ QYRALTVP EL TQQMWDSKNM MCAADPRHGR YLTASAMFRG KMSTKEVDEQ MINVQKNSS YFVEWIPNNV KSSVCDIAPR GLSMAS TFIG NSTSIQEMFR RVSEQFTAMF RRKAFLHWYT GEGMDEMEFT EAESNMNDLV SEYQQYQDAT ADEEGEYEDE EEEEPEHG YE
Source	Mammalian cell
Target Names	TUBB3
Protein Names	Recommended name: Tubulin beta-3 chain Alternative name(s): Beta-3-tubulin
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.