



Recombinant *Caenorhabditis elegans* Protein Wnt-2 (wnt-2)

Product Code	CSB-EP026133CXY-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34889
Product Type	Recombinant Protein
Immunogen Species	<i>Caenorhabditis elegans</i>
Purity	≥85% (SDS-PAGE)
Sequence	AGH NVKPVCELPG LSPGQAQVCE LFKDHMPAVS IGAQNAIQEC QRQFTGHRWN CSTHYSTGML GPIHKMATRE AAFYAILS A GVTHEIGRRC KQGLLTSCGC SDETKPKNVP TDWSWGGCGD NVEYGYKFSR DFIDIREKEH DPKRNHDNGR SLMNRRNNEA GRKILKRHRK PKCKCHGVSG ACNMKTCWMQ LPSMEQVGKI LRNKYDKAIR VQINDRGNLQ LLADEATKER KTRALPTDLV FMDDSPDYCR FDRHSGTLGT EGRVCKRGS G GAEGCDSLCC GRGYNTYTQE VKSKCNCKFE WCKVVCQTC NNVTQVDICK
Source	E.coli
Target Names	cwn-2
Protein Names	Recommended name: Protein Wnt-2
Expression Region	38-360
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.