



# Recombinant Pyruvate dehydrogenase E1 component subunit beta (pdhB)

<b>Product Code</b>	CSB-MP017717ABI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P35488
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Acholeplasma laidlawii
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAIITLLEAI NQAIDQAMEK DESIVVFGED AGFEGGVFRV TAGLQKKYGE TRVFDTPIAE SAIVGSVGM AINGLKPIAE IQFDGFIFPG YTDLVTHAAR MRNRSRGQFT VPMVLRLPHG GGIRALEHHS EALEVLFGSI PGLKVVTPST PYDAKGLLLA AINDPDPVVF LEPKRIYRAG KQEVPAEMYE IPIGKAKVVK QGDTMTVVAW GSIVREVEKA VKLVEAEGIS VEIIDLRTIS PIDEETILNS VKKTGKFMVV TEAVKSYGPA AELITMVNEK AFFHLEAAPV RFTGFDITVP LARGEHYHFP QPEKIAAYIR KLAKARP
<b>Source</b>	Mammalian cell
<b>Target Names</b>	pdhB
<b>Protein Names</b>	Recommended name: Pyruvate dehydrogenase E1 component subunit beta EC= 1.2.4.1
<b>Expression Region</b>	1-327
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.