



Recombinant Giardia intestinalis Triosephosphate isomerase

Product Code	CSB-EP024102GFW-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P36186
Product Type	Recombinant Protein
Immunogen Species	Giardia intestinalis (Giardia lamblia)
Purity	>85% (SDS-PAGE)
Sequence	MPARRPFIGG NFKCNGSLDF IKSHVAAIAA HKIPDSVDVV IAPSAVHLST AIAANTSKQL RIAAQNVYLE GNGAWTGETS VEMLQDMGLK HVIVGHSERR RIMGETDEQS AKKAKRALEK GMTVIFCVGE TLDERKANRT MEVNIAQLEA LGKELGESKM LWKEVVIAYE PVWSIGTGVV ATPEQAEVH VGLRKWFAEK VCAEGAQHIR IYGGSSANGS NCEKLGQCPN IDGFLVGGAS LKPEFMTMID ILTKTRT
Source	E.coli
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	1-257
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.