



# Recombinant Mycoplasma hominis Elongation factor Tu (tuf)

<b>Product Code</b>	CSB-EP025342MWJ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P22679
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mycoplasma hominis (strain ATCC 23114 / NBRC 14850 / NCTC 10111 / PG21)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MAKLDFDRSK PHVNIGTIGH VDHGKTTTLTA AIATVLAKKG LAEARDYASI DNAPEEKARG ITINTSHIEY QTEKRHYAHV DCPGHADYVK NMITGAAQMD GAILVVAATD GPMPQTREHI LLARQVGVPK IVVFLNKIDM FKDDEREEMV GLVEMDVRS LSEYGFDDGN APIIAGSALK ALQGDPEYEK GILELMDAVD TYIEEPKRET DKPFLMAVED VFTITGRGTV ATGRVERGVL QLNEEVEIVG LKPTKKT VVT GIEMFRKNLK EAQAGDNAGL LLRGIDRSEV ERGQVLAKPK TIVPHTQFEA TVYVLKKEEG GRHTPPFFHNY KPQFYFRTTD VTGGIEFKPG REMVVPGDNV ELTVTLI API AIEEGTKFSI REGGRTVGAG SVTKILK
<b>Source</b>	E.coli
<b>Target Names</b>	tuf
<b>Protein Names</b>	Recommended name: Elongation factor Tu Short name= EF-Tu
<b>Expression Region</b>	1-397
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.