



Recombinant Mycoplasma hominis Elongation factor Tu (tuf)

Product Code	CSB-MP025342MWJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P22679
Product Type	Recombinant Protein
Immunogen Species	Mycoplasma hominis (strain ATCC 23114 / NBRC 14850 / NCTC 10111 / PG21)
Purity	>85% (SDS-PAGE)
Sequence	MAKLDFDRSK PHVNIGTIGH VDHGKTTTLTA AIATVLAKKG LAEARDYASI DNAPEEKARG ITINTSHIEY QTEKRHYAHV DCPGHADYVK NMITGAAQMD GAILVVAATD GPMPQTREHI LLARQVGVPK IVVFLNKIDM FKDDEREEMV GLVEMDVRS LSEYGFDDGN APIIAGSALK ALQGDPEYEK GILELMDAVD TYIEEPKRET DKPFLMAVED VFTITGRGTV ATGRVERGVL QLNEEVEIVG LKPTKKT VVT GIEMFRKNLK EAQAGDNAGL LLRGIDRSEV ERGQVLAKPK TIVPHTQFEA TVYVLKKEEG GRHTPFFHNY KPQFYFRTTD VTGGIEFKPG REMVVPGDNV ELTVTLI API AIEEGTKFSI REGGRTVGAG SVTKILK
Source	Mammalian cell
Target Names	tuf
Protein Names	Recommended name: Elongation factor Tu Short name= EF-Tu
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.