



Recombinant *Drosophila melanogaster* Tubulin gamma-1 chain (gammaTub23C)

Product Code	CSB-BP025334DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P23257
Product Type	Recombinant Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MPSEIITLQL GQCGNQIGFE FWKRLCLEHG ISPSGVLEDF ANDGLDRKDV FFYQADDDHY IPRAVLLDLE PRVINTIMGS VYSKLYNPEN VYLSKHGGGA GNNWASGYSQ GEKLQEEVFD IIDREADGSD SLEGFILCHS IAGGTGSGMG SFIMERLADR YPKKLIQTFS VFPNQDEISD VVVQPYNMML TLKRLTTAAD SVVVDNTAL NRIACDRLHI QNPSFSQINN LVSTIMSVST TTLRYPYMM NNLIGLTAPL IPTPQLHFLM TGYTPLTSDS DIHTQQLVNV RKTTVLDVMR RLLQPKNMMV STGPDKSNHH CYISILNIIQ GEVDPTQVHK SLQRIRDRKM AQFIPWGPTS IQVALSRSSP YVQSNHRVSG LMLANHTSIC SLFERALNQY DKLRKRGAFD DQFRREDIFK DDLNELDESR ETVDCLVQEQ EAATREYMQ FSVKGNGPV DSKSEDSRSV TSAGS
Source	Baculovirus
Target Names	gammaTub23C
Protein Names	Recommended name: Tubulin gamma-1 chain Alternative name(s): Gamma-1-tubulin
Expression Region	1-475
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.