



Recombinant Solanum tuberosum TATA-box-binding protein (TBP)

Product Code	CSB-EP023239FIG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P26357
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MADQGLEGSQ PVDLTKHPSG IVPTLQNIIVS TVNLDCKLDL KAIALQARNA EYNPKRFAAV IMRIREPKTT ALIFASGKMV CTGAKSEQQS KLAARKYARI IQKLGFPKAF KDFKIQNIVG SCDVKFPIRL EGLAYAHGAF SSYEPELFPG LIYRMKQPKI VLLIFVSGKI VITGAKVRDE TYTAFENIYP VLTEFRKNQQ
Source	E.coli
Target Names	TBP
Protein Names	Recommended name: TATA-box-binding protein Alternative name(s): TATA sequence-binding protein Short name= TBP TATA-binding factor TATA-box factor Transcription initiation factor TFIID TBP subunit
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.